

Product Information Sheet:

Animal Product-Free LB Broth (Miller), (Lennox), (Luria)



AssayCell Technologies
36 Chesnut St.
Avenel, NJ 07001

Ordering: 1-732-595-7680 or sales@assaycell.com

Product Information Sheet: *Private Labeling Users Only*

Product Name:

APF LB Broth (Miller)
APF LB Broth (Lennox)
APF LB Broth (Luria)

Product Description:

The APF LB Media are animal-product-free versions of the industry standard LB Broth media. These media are prepared using Atholate™, AthenaES' proprietary plant protein hydrolysate. (See *Technical Brief: APF LB Broth.*)

APF LB Media are certified animal-product-free.

Appearance: tan, free-flowing homogenous powder.

Intended Use:

APF LB Media are used for the cultivation of *Escherichia coli*, *Salmonella*, *Pseudomonas*, *Bacillus* and other eubacteria species.

Catalog Number: 0131 (Luria)
0132 (Lennox)
0133 (Miller)

Formulations:

Compositions of APF LB Broths (g/L)			
Ingredient	Miller	Lennox	Luria
Yeast Extract	5	5	5
NaCl	10	5	0.5
Atholate™	15	15	15

Bulk packaging is available. Please inquire.

Availability and Pricing (in US dollars):

For distribution under the Athena Enzyme Systems™ brand.
See *Product Pricing Schedule* below.

Safety: Full MSDS available upon request.

Handling and Storage:

Handling: Information for Safe Handling:
No special measures required.
Information about protection against explosions and fires:
No special measures required.

Storage: Requirements to be met by storerooms and receptacles:
No special requirements.
Information about storage in one common storage facility:
Store away from oxidizing agents.
Further information about storage conditions: None.

Hazard Identification:

Hazard description:
This product contains no hazardous constituents, or the concentration of all chemical constituents are below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC, and 93/112/EC.

NEPA ratings (scale 0-4):
Health = 0
Fire = 0
Reactivity = 0

Atholate™ is a trademark of Athena Environmental Sciences, Inc.

Product Pricing Schedule

Product Code	Product Name	Unit Size	Suggested Retail Price	Discount Level 1	Discount Level 2	Discount Level 3	Discount Level 4	Bulk (25 kg)
Expression Media Products			1 Unit	<50 Units	51-100 Units	101-150 Units	>150 Units	
0131	APF LF Luria	500 g	\$ 40.00	\$ 26.00	\$ 24.70	\$ 23.40	\$ 22.10	\$ 1,000.00
0132	APF LB Lennox	500 g	\$ 37.00	\$ 24.05	\$ 22.85	\$ 21.65	\$ 20.44	\$ 925.00
0133	APF LB Miller	500 g	\$ 35.00	\$ 22.75	\$ 21.61	\$ 20.48	\$ 19.34	\$ 875.00
0131-S	APF LF Luria, 1L Singles	Box 10	\$ 40.00	\$ 26.00	\$ 24.70	\$ 23.40	\$ 22.10	
0132-S	APF LB Lennox, 1L Singles	Box 10	\$ 37.00	\$ 24.05	\$ 22.85	\$ 21.65	\$ 20.44	
0133-S	APF LB Miller, 1L Singles	Box 10	\$ 35.00	\$ 22.75	\$ 21.61	\$ 20.48	\$ 19.34	