

Product Information Sheet:

Atholate™ (Animal Product-Free **casein hydrolysate replacement**)



AssayCell Technologies
36 Chesnut St
Avenel, NJ 07001

Ordering: 1-732-595-7680 or sales@assaycell.com

Product Information Sheet: *Private Labeling Users Only*

Product Name: Atholate™

Product Description:

Atholate™ is a proprietary blend of plant protein hydrolysates designed to replace casein-derived hydrolysates by matching the nutritional and performance characteristics of casein hydrolysates. (See *Technical Brief: APF LB Broth*.)

Atholate™ is certified animal-product-free.

Appearance: tan, free-flowing homogenous powder.

Intended Use:

Atholate™ is used to prepare bacteriological media for the cultivation of a wide range of microbial species. Use 1.5g Atholate™ per 1.0g of casein hydrolysate.

Catalog Number:

0134-010
0134-025
0134-050

Availability and Pricing (in US dollars):

For private labeling or formulation in generic media products. Product labels must list Atholate™ as an ingredient. Pricing is based on amount of product ordered with an institutional purchase order. New customers must prepay on the first order if the purchase order is for more than \$1,000 distribution under the Athena Enzyme Systems™ brand.

See *Product Pricing Schedule* below.

Package sizes greater than 10kg are available. Please inquire.

Safety: Full MSDS available upon request.

Handling and Storage:

Handling: *Information for Safe Handling:*

No special measures required.

Information about protection against explosions and fires:

No special measures required.

Storage: *Requirements to be met by storerooms and receptacles:*

No special requirements.

Information about storage in one common storage facility:

Store away from oxidizing agents.

Further information about storage conditions: None.

Hazard Identification:

Hazard description:

This product contains no hazardous constituents, or the concentration of all chemical constituents are below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC, and 93/112/EC.

NFPA ratings (scale 0-4):

Health = 0

Fire = 0

Reactivity = 0

Atholate™ is a trademark of Athena Environmental Sciences, Inc.

Product Pricing Schedule

Package/Order Sizes	Price per kg			
	<500 kg	500-1000 kg	1000-2000 kg	>2000 kg
10 kg	\$ 43.00	\$40.85	\$39.78	\$38.70
25 kg	\$ 42.14	\$40.03	\$38.98	\$37.93
50 kg	\$41.28	\$39.22	\$38.18	\$37.15
500 g Bottle	Price per 500 g Container			
	<50	51-100	101-150	>150
500 g Bottle	\$42.25	\$40.14	\$38.03	\$35.91