

2018 Winter Special

BeadBlaster™ 24 Place Homogenizer



Offer Valid Through: 06.30.18

Promo Code: WP2400

BEADBLASTER™ 24

MICROTUBE HOMOGENIZER

\$6495*
Valid thru. 6/30/18

~~\$7995~~

- + Extremely powerful mixing action, for **lysis, grinding or homogenization**
- + Sealed 2.0ml tubes - no cross-contamination
- + Faster & more effective than tissue grinders
- + Ideal for extraction of DNA, RNA and Proteins from all sample types



Technical Data:

Capacity:	24 x 2.0ml tubes
Speed:	4.00 to 7.00 m/s, 0.05 increments
Cycle time:	1 - 90 sec, 1 sec increments
Pause between cycles:	0-2 min, 1 sec increments
Max. cycles per program:	10
Program memory:	up to 50
Accel. / Decel. time:	<2 sec / <2 sec
Noise level:	<68db
Dimensions (WxDxH):	11 x 14.2 x 15.2 in / 28 x36 x38.5 cm
Weight:	55 lb / 25kg
Electrical:	100V to 240V, 50-60Hz, 600W
Warranty:	2 years

New! BeadBug™ 6



BeadBug™ 6 Microtube Homogenizer for 6x2.0ml tubes or 2x5.0ml tubes* (with optional adapters)

Zirconium



Triple-Pure™: Acid washed, heat treated & tested free of nuclease

Triple-Pure™ (Molecular Biology Grade) Zirconium Beads

Ordering Information:

Item No.	Description:	Price:
SKU:3366	BeadBlaster™ 24 Microtube Homogenizer* with 24 x 2.0ml tube holder	\$6495
SKU:3368	Optional tube holder, 12 x 5ml	\$345.00
Prefilled Tube Kits:	Include 2.0ml non-skirted screw cap tubes with caps, sealing rings and beads as specified below:	
SKU:1613	Standard glass beads, Ø 0.1mm, 50pk.	\$104.50
SKU:1614	Standard glass beads, Ø 0.5mm, 50pk.	\$104.50
SKU:1615	Standard glass beads, Ø 1.0mm, 50pk.	\$104.50
SKU:1616	Triple-Pure, high impact zirconium beads, Ø 0.1mm, 50pk.	\$146.00
SKU: 1617	Triple-Pure, high impact zirconium beads, Ø 0.5mm, 50pk.	\$146.00
	Triple-Pure, high impact zirconium beads, Ø 1.0mm, 50pk.	\$146.00
SKU:1619	Triple-Pure, high impact zirconium beads, Ø 1.5mm, 50pk.	\$146.00
SKU:1620	Triple-Pure, high impact zirconium beads, Ø 3.0mm, 50pk.	\$146.00
SKU:1621	Starter Kit, 10 ea. of Ø 0.1, 0.5, 1.0, 1.5 & 3.0 Triple-Pure, high impact zirconium	\$146.00
SKU:1625	One ceria zirconium oxide satellite, Ø 6.0mm, 50pk.	\$146.00
SKU:1622	Stainless steel beads, Ø 2.8mm, 50pk.	\$146.00
SKU:1623	Garnet shards with 1 zirconium bead Ø 6.0mm, 50pk.	\$146.00
SKU:1624	Mixed format, 0.1mm silica, 1.4mm zirconium, 4mm silica, 50pk	\$146.00
SKU:1599	Empty tubes (2.0ml) pack of 1000	\$154.00
SKU:1601	Reinforced empty tubes (2.0ml) pack of 500	\$189.00
SKU: 1626	5ml tubes (snap cap) with stainless beads, Ø 2.8mm, 50pk	\$195.00

Additional bead types are also available (in both bulk and prefilled kits). Please contact Benchmark for a complete listing.



Video Demonstration